

For Research Use Only

MMP8 Recombinant antibody, PBS Only (Capture/Detector)



Catalog Number: 82885-5-PBS

Basic Information

Catalog Number:
82885-5-PBS

Size:
1mg/ml

Source:
Rabbit

Isotype:
IgG

GenBank Accession Number:
BC074988

GeneID (NCBI):
4317

UNIPROT ID:
P22894

Full Name:
matrix metalloproteinase 8
(neutrophil collagenase)

Calculated MW:
467 aa, 53 kDa

Purification Method:
Protein A purification

CloneNo.:
230030H8

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP00013-1, MMP8 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82885-5-PBS. Detection antibody: 82885-4-PBS. Standard: Eg0378. Range: 0.156-10 ng/mL



Sandwich ELISA standard curve of MP00013-2, MMP8 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82885-3-PBS. Detection antibody: 82885-5-PBS. Standard: Eg0378. Range: 0.156-10 ng/mL