

For Research Use Only

IFT88 Recombinant antibody

Catalog Number: 82886-2-RR



Basic Information

Catalog Number: 82886-2-RR	GenBank Accession Number: BC030776	Purification Method: Protein A purification
Size: 1000 µg/ml	GeneID (NCBI): 8100	CloneNo.: 230074E12
Source: Rabbit	UNIPROT ID: Q13099	Recommended Dilutions: WB 1:5000-1:50000 IF 1:200-1:800
Isotype: IgG	Full Name: intraflagellar transport 88 homolog (Chlamydomonas)	
Immunogen Catalog Number: AG4980	Calculated MW: 94 kDa	
	Observed MW: 94 kDa	

Applications

Tested Applications: WB, IF/ICC, ELISA	Positive Controls: WB : HEK-293 cells, A549 cells, PC-3 cells IF : Starvation treated hTERT-RPE1 cells,
Species Specificity: Human	

Background Information

Intraflagellar transport (IFT), mediated by molecular motors and IFT particles, is an important transport process that occurs in the cilium and has been shown to be essential for the assembly and maintenance of cilia and flagella in many organisms. IFT88 (intraflagellar transport protein 88; also known as TG737 or TTC10) is a component of IFT particles and required for cilium biogenesis. Defects in IFT88/Tg737 lead to polycystic kidney disease (11062270). IFT88 localizes to spindle poles during mitosis and is required for spindle orientation in mitosis (21441926).

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 82886-2-RR (IFT88 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.



Immunofluorescent analysis of (4% PFA) fixed Starvation treated hTERT-RPE1 cells using IFT88 antibody (82886-2-RR, Clone: 230074E12) at dilution of 1:400 and CoraLite®488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) (SA00013-2), Acetyl-Tubulin (Lys40) antibody (66200-1-Ig, Clone: 7E5H8, red).