

For Research Use Only

# TRIM62 Recombinant antibody, PBS Only (Capture)

Catalog Number:82892-3-PBS



## Basic Information

Catalog Number: 82892-3-PBS	GenBank Accession Number: BC001222	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 55223	CloneNo.: 230195C11
Source: Rabbit	UNIPROT ID: Q9BVG3	
Isotype: IgG	Full Name: tripartite motif-containing 62	
Immunogen Catalog Number: AG7143	Calculated MW: 54 kDa	

## Applications

Tested Applications:  
Species Specificity:  
human

## Background Information

### Storage

Storage:  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**  
Storage Buffer:  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP01312-2, TRIM62 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82892-3-PBS. Detection antibody: 82892-4-PBS. Standard: Ag7143. Range: 1.25-20 ng/mL