

For Research Use Only

# TRIM62 Recombinant antibody, PBS Only (Detector)

Catalog Number:82892-4-PBS



## Basic Information

Catalog Number:

82892-4-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG7143

GenBank Accession Number:

BC001222

GeneID (NCBI):

55223

UNIPROT ID:

Q9BVG3

Full Name:

tripartite motif-containing 62

Calculated MW:

54 kDa

Purification Method:

Protein A purification

CloneNo.:

230195C10

## Applications

Tested Applications:

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP01312-2, TRIM62 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82892-3-PBS. Detection antibody: 82892-4-PBS. Standard: Ag7143. Range: 1.25-20 ng/mL