

For Research Use Only

NSUN2 Recombinant antibody, PBS Only (Capture)



Catalog Number: 82894-6-PBS

Basic Information

Catalog Number: 82894-6-PBS	GenBank Accession Number: BC001041	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 54888	CloneNo.: 230186F11
Source: Rabbit	UNIPROT ID: Q08J23	
Isotype: IgG	Full Name: NOL1/NOP2/Sun domain family, member 2	
Immunogen Catalog Number: AG14791	Calculated MW: 767 aa, 86 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00112-2, NSUN2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82894-6-PBS. Detection antibody: 82894-5-PBS. Standard: Ag14791. Range: 0.625-20 ng/mL