

For Research Use Only

CHD5 Recombinant antibody, PBS Only



Catalog Number: 82896-1-PBS

Basic Information

Catalog Number:

82896-1-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG18621

GenBank Accession Number:

BC148803

GeneID (NCBI):

26038

UNIPROT ID:

Q8TD10

Full Name:

chromodomain helicase DNA binding protein 5

Calculated MW:

1954 aa, 223 kDa

Purification Method:

Protein A purification

CloneNo.:

230093A6

Applications

Tested Applications:

Indirect ELISA

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data