

For Research Use Only

CHD5 Recombinant antibody

Catalog Number: 82896-1-RR



Basic Information

Catalog Number: 82896-1-RR	GenBank Accession Number: BC148803	Purification Method: Protein A purification
Size: 1000 µg/ml	GeneID (NCBI): 26038	CloneNo.: 230093A6
Source: Rabbit	UNIPROT ID: Q8TD10	
Isotype: IgG	Full Name: chromodomain helicase DNA binding protein 5	
Immunogen Catalog Number: AG18621	Calculated MW: 1954 aa, 223 kDa	

Applications

Tested Applications:
ELISA

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.

Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data