

For Research Use Only

FANCM Recombinant antibody, PBS Only (Capture)



Catalog Number: 82897-3-PBS

Basic Information

Catalog Number: 82897-3-PBS	GenBank Accession Number: BC140776	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 57697	CloneNo.: 230233A10
Source: Rabbit	UNIPROT ID: Q8IYD8	
Isotype: IgG	Full Name: Fanconi anemia, complementation group M	
Immunogen Catalog Number: AG34447		

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00057-3, FANCM Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82897-3-PBS. Detection antibody: 82897-4-PBS. Standard: Ag34447. Range: 0.313-20 ng/mL.