

For Research Use Only

# FANCM Recombinant antibody, PBS Only (Detector)



Catalog Number: 82897-4-PBS

## Basic Information

Catalog Number:

82897-4-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG34447

GenBank Accession Number:

BC140776

GeneID (NCBI):

57697

UNIPROT ID:

Q8IYD8

Full Name:

Fanconi anemia, complementation

group M

Purification Method:

Protein A purification

CloneNo.:

230233B4

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00057-2, FANCM Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82897-6-PBS. Detection antibody: 82897-4-PBS. Standard: Ag34447. Range: 0.313-20 ng/mL



Cytometric bead array standard curve of MP00057-3, FANCM Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82897-3-PBS. Detection antibody: 82897-4-PBS. Standard: Ag34447. Range: 0.313-20 ng/mL