

For Research Use Only

# FANCM Recombinant antibody, PBS Only (Capture)



Catalog Number: 82897-6-PBS

## Basic Information

|   |  |   |
|---|--|---|
| <b>Catalog Number:</b><br>82897-6-PBS       | <b>GenBank Accession Number:</b><br>BC140776                 | <b>Purification Method:</b><br>Protein A purification |
| <b>Size:</b><br>1mg/ml                      | <b>GeneID (NCBI):</b><br>57697                               | <b>CloneNo.:</b><br>230233G9                          |
| <b>Source:</b><br>Rabbit                    | <b>UNIPROT ID:</b><br>Q8IYD8                                 |   |
| <b>Isotype:</b><br>IgG                      | <b>Full Name:</b><br>Fanconi anemia, complementation group M |   |
| <b>Immunogen Catalog Number:</b><br>AG34447 |  |   |

## Applications

**Tested Applications:**  
Indirect ELISA, Cytometric bead array

**Species Specificity:**  
Human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00057-1, FANCM Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82897-6-PBS. Detection antibody: 82897-5-PBS. Standard: Ag34447. Range: 0.313-20 ng/mL



Cytometric bead array standard curve of MP00057-2, FANCM Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82897-6-PBS. Detection antibody: 82897-4-PBS. Standard: Ag34447. Range: 0.313-20 ng/mL