

For Research Use Only

HLA-DMA Recombinant antibody, PBS Only (Capture)



Catalog Number: 82899-3-PBS

Basic Information

Catalog Number:

82899-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG2452

GenBank Accession Number:

BC026279

GeneID (NCBI):

3108

UNIPROT ID:

Q31604

Full Name:

major histocompatibility complex,
class II, DM alpha

Calculated MW:

261 aa, 29 kDa

Purification Method:

Protein A purification

CloneNo.:

230155C7

Applications

Tested Applications:

Cytometric bead array, ELISA

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Standard curve of MP00018-3 measured by cytometric bead array. Capture antibody: 82899-3-PBS. Detection antibody: 82899-1-PBS. Range: 0.625-20 ng/mL