

For Research Use Only

# Synaptophysin; SYP Recombinant antibody, PBS Only (Capture)



Catalog Number: 82900-5-PBS

## Basic Information

Catalog Number: 82900-5-PBS	GenBank Accession Number: BC064550	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 6855	CloneNo.: 240162D5
Source: Rabbit	UNIPROT ID: P08247	
Isotype: IgG	Full Name: synaptophysin	
Immunogen Catalog Number: AG11781	Calculated MW: 313 aa, 34 kDa	

## Applications

**Tested Applications:**  
Indirect ELISA, Cytometric bead array

**Species Specificity:**  
Human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00345-2, Synaptophysin; SYP Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82900-5-PBS. Detection antibody: 82900-4-PBS. Standard: Ag11781. Range: 0.625-80 ng/mL.