

For Research Use Only

# BAZ2A Recombinant antibody, PBS Only (Detector)

Catalog Number: 82910-3-PBS



## Basic Information

Catalog Number:

82910-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG33963

GenBank Accession Number:

NM\_013449.4

GeneID (NCBI):

11176

UNIPROT ID:

Q9UIF9

Full Name:

bromodomain adjacent to zinc finger domain, 2A

Calculated MW:

211 kDa

Purification Method:

Protein A purification

CloneNo.:

230205F10

## Applications

Tested Applications:

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP01313-2, BAZ2A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82910-4-PBS. Detection antibody: 82910-3-PBS. Standard: Ag33963. Range: 1.25-20 ng/mL