

For Research Use Only

TBCD Recombinant antibody, PBS Only



Catalog Number: 82911-2-PBS

Basic Information

Catalog Number: 82911-2-PBS	GenBank Accession Number: BC003094	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 6904	CloneNo.: 230209D7
Source: Rabbit	UNIPROT ID: Q9BTW9	
Isotype: IgG	Full Name: tubulin folding cofactor D	
Immunogen Catalog Number: AG6310	Calculated MW: 133 kDa	

Applications

Tested Applications:
Indirect ELISA
Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Biolayer interferometry (BLI) kinetic assays of 82911-2-RR against Human TBCD were performed. The affinity constant is 2.07 nM.