

For Research Use Only

ZC3H11A Recombinant antibody, PBS Only



Catalog Number: 82912-5-PBS

Basic Information

Catalog Number: 82912-5-PBS	GenBank Accession Number: BC014268	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 9877	CloneNo.: 230191D2
Source: Rabbit	UNIPROT ID: O75152	
Isotype: IgG	Full Name: zinc finger CCCH-type containing 11A	
Immunogen Catalog Number: AG23424		

Applications

Tested Applications:
IF, Indirect ELISA
Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed HepG2 cells using ZC3H11A antibody (82912-5-RR, Clone: 230191D2) at dilution of 1:400 and CoraLite@488-Conjugated AffiniPure Goat Anti-Mouse IgG(H+L) (SA00013-1), CL594-Phalloidin (red). This data was developed using the same antibody clone with 82912-5-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 82912-5-RR against Human ZC3H11A were performed. The affinity constant is 0.328 nM.