

For Research Use Only

# MYOF Recombinant antibody

Catalog Number: 82917-7-RR



## Basic Information

Catalog Number:

82917-7-RR

Size:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_133337

GeneID (NCBI):

26509

UNIPROT ID:

Q9NZM1

Full Name:

myoferlin

Calculated MW:

235 kDa

Purification Method:

Protein A purification

CloneNo.:

230213B6

Recommended Dilutions:

IF/ICC 1:125-1:500

## Applications

Tested Applications:

IF/ICC, ELISA

Species Specificity:

Human

Positive Controls:

IF/ICC : HeLa cells,

## Background Information

### Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed HeLa cells using MYOF antibody (82917-7-RR, Clone: 230213B6 ) at dilution of 1:250 and CoraLite®488-Conjugated AffiniPure Goat Anti-Mouse IgG(H+L) (SA00013-1).