

For Research Use Only

PICK1 Recombinant antibody, PBS Only (Capture)



Catalog Number: 82919-3-PBS

Basic Information

| | | |
|--|---|---|
| Catalog Number: 82919-3-PBS | GenBank Accession Number: BC017561 | Purification Method: Protein A purification |
| Size: 1 mg/ml | GeneID (NCBI): 9463 | CloneNo.: 230185A9 |
| Source: Rabbit | UNIPROT ID: Q9NRD5 | |
| Isotype: IgG | Full Name: protein interacting with PRKCA 1 | |
| Immunogen Catalog Number: AG1443 | Calculated MW: 47 kDa | |

Applications

Tested Applications:
Cytometric bead array, ELISA

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Standard curve of MP00068-2 measured by cytometric bead array. Capture antibody: 82919-3-PBS. Detection antibody: 82919-1-PBS. Range: 0.625-20 ng/mL