

For Research Use Only

NFX1 Recombinant antibody

Catalog Number:82924-1-RR



Basic Information

Catalog Number:

82924-1-RR

Size:

1000 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG9448

GenBank Accession Number:

BC012151

GeneID (NCBI):

4799

UNIPROT ID:

Q12986

Full Name:

nuclear transcription factor, X-box binding 1

Calculated MW:

1120 aa, 124 kDa

Purification Method:

Protein A purification

CloneNo.:

230211D9

Recommended Dilutions:

IF/ICC 1:200-1:800

Applications

Tested Applications:

IF/ICC, FC (Intra), ELISA

Species Specificity:

human

Positive Controls:

IF/ICC : HeLa cells,

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:



T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1×10^6 PC-3 cells were intracellularly stained with 0.25 μ g NFX1 Recombinant antibody (82924-1-RR, Clone:230211D9) and APC-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)(red), or 0.25 μ g Isotype Control (blue). Cells were fixed and permeabilized with True-Nuclear Transcription Factor Buffer Set.

Immunofluorescent analysis of (4% PFA) fixed HeLa cells using NFX1 antibody (82924-1-RR, Clone: 230211D9) at dilution of 1:400 and CoraLite® 488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) (SA00013-2).