

For Research Use Only

NFX1 Recombinant antibody, PBS Only



Catalog Number: 82924-2-PBS

Basic Information

Catalog Number: 82924-2-PBS	GenBank Accession Number: BC012151	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 4799	CloneNo.: 230211H4
Source: Rabbit	UNIPROT ID: Q12986	
Isotype: IgG	Full Name: nuclear transcription factor, X-box binding 1	
Immunogen Catalog Number: AG9448	Calculated MW: 1120 aa, 124 kDa	
	Observed MW: 124 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 82924-2-RR (NFX1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 82924-2-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 82924-2-RR against Human NFX1 were performed. The affinity constant is below 1 pM.