

For Research Use Only

NFX1 Recombinant antibody, PBS Only (Detector)



Catalog Number:82924-6-PBS

Basic Information

Catalog Number: 82924-6-PBS	GenBank Accession Number: BC012151	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 4799	CloneNo.: 230211H6
Source: Rabbit	UNIPROT ID: Q12986	
Isotype: IgG	Full Name: nuclear transcription factor, X-box binding 1	
Immunogen Catalog Number: AG9448	Calculated MW: 1120 aa, 124 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00127-2, NFX1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82924-3-PBS. Detection antibody: 82924-6-PBS. Standard: Ag9448. Range: 0.313-20 ng/mL



Cytometric bead array standard curve of MP00127-3, NFX1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82924-5-PBS. Detection antibody: 82924-6-PBS. Standard: Ag9448. Range: 0.625-20 ng/mL