

For Research Use Only

MAPKAPK2 Recombinant antibody

Catalog Number: 82929-2-RR



Basic Information

Catalog Number: 82929-2-RR	GenBank Accession Number: BC036060	Purification Method: Protein A purification
Size: 700 µg/ml	GeneID (NCBI): 9261	CloneNo.: 230197D6
Source: Rabbit	UNIPROT ID: P49137	Recommended Dilutions: WB 1:2000-1:14000
Isotype: IgG	Full Name: mitogen-activated protein kinase- activated protein kinase 2	
Immunogen Catalog Number: AG34196	Calculated MW: 400 aa, 46 kDa	
	Observed MW: 42-46 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : A549 cells, HeLa cells, HEK-293T cells
Species Specificity: Human	

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 82929-2-RR (MAPKAPK2 antibody) at dilution of 1:7000 incubated at room temperature for 1.5 hours.