

For Research Use Only

RORA Recombinant antibody, PBS Only (Capture)



Catalog Number: 82930-6-PBS

Basic Information

Catalog Number:

82930-6-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0944

GenBank Accession Number:

BC008831

GeneID (NCBI):

6095

UNIPROT ID:

P35398

Full Name:

RAR-related orphan receptor A

Calculated MW:

59 kDa

Purification Method:

Protein A purification

CloneNo.:

230217E4

Applications

Tested Applications:

Cytometric bead array, ELISA

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Standard curve of MP00055-2 measured by cytometric bead array. Capture antibody: 82930-6-PBS. Detection antibody: 82930-2-PBS. Range: 0.625-20 ng/mL



Standard curve of MP00055-3 measured by cytometric bead array. Capture antibody: 82930-6-PBS. Detection antibody: 82930-4-PBS. Range: 0.625-20 ng/mL