

For Research Use Only

EIF2S1 Recombinant antibody, PBS Only (Detector)



Catalog Number: 82936-2-PBS

Basic Information

Catalog Number: 82936-2-PBS	GenBank Accession Number: BC002513	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 1965	CloneNo.: 230245A2
Source: Rabbit	UNIPROT ID: P05198	
Isotype: IgG	Full Name: eukaryotic translation initiation factor 2, subunit 1 alpha, 35kDa	
Immunogen Catalog Number: AG1645	Calculated MW: 36 kDa	

Applications

Tested Applications:
Cytometric bead array, ELISA

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Standard curve of MP00076-1 measured by cytometric bead array. Capture antibody: 82936-3-PBS. Detection antibody: 82936-2-PBS. Range: 0.313-20 ng/mL