

For Research Use Only

# EIF2S1 Recombinant antibody, PBS Only (Capture)



Catalog Number: 82936-3-PBS

## Basic Information

Catalog Number:

82936-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG1645

GenBank Accession Number:

BC002513

GeneID (NCBI):

1965

UNIPROT ID:

P05198

Full Name:

eukaryotic translation initiation factor 2, subunit 1 alpha, 35kDa

Calculated MW:

36 kDa

Purification Method:

Protein A purification

CloneNo.:

230245B4

## Applications

Tested Applications:

Cytometric bead array, ELISA

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Standard curve of MP00076-1 measured by cytometric bead array. Capture antibody: 82936-3-PBS. Detection antibody: 82936-2-PBS. Range: 0.313-20 ng/mL