

For Research Use Only

EIF2S1 Recombinant antibody, PBS Only (Detector)



Catalog Number: 82936-7-PBS

Basic Information

Catalog Number:

82936-7-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG1645

GenBank Accession Number:

BC002513

GeneID (NCBI):

1965

UNIPROT ID:

P05198

Full Name:

eukaryotic translation initiation factor 2, subunit 1 alpha, 35kDa

Calculated MW:

36 kDa

Purification Method:

Protein A purification

CloneNo.:

230245H5

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00076-3, EIF2S1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82936-6-PBS. Detection antibody: 82936-7-PBS. Standard: Ag1645. Range: 1.25-20 ng/mL.