

For Research Use Only

CAMP Recombinant antibody, PBS Only (Capture)



Catalog Number: 82940-4-PBS

Basic Information

Catalog Number: 82940-4-PBS	GenBank Accession Number: BC055089	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 820	CloneNo.: 230157E9
Source: Rabbit	UNIPROT ID: P49913	
Isotype: IgG	Full Name: cathelicidin antimicrobial peptide	
Immunogen Catalog Number: AG2622	Calculated MW: 170 aa, 19 kDa	

Applications

Tested Applications:
Cytometric bead array, ELISA

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Standard curve of MP00042-2 measured by sandwich ELISA. Capture antibody: 82940-4-PBS. Detection antibody: 82940-6-PBS. Range: 0.39-25 ng/mL