

For Research Use Only

# PDGFRB Recombinant antibody, PBS Only (Detector)



Catalog Number: 82943-4-PBS

## Basic Information

Catalog Number:

82943-4-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC032224

GeneID (NCBI):

5159

Full Name:

platelet-derived growth factor receptor, beta polypeptide

Calculated MW:

1106 aa, 124 kDa

Purification Method:

Protein A purification

CloneNo.:

230214E10

## Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Cytometric bead array standard curve of MP00006-1, PDGFRB Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82943-2-PBS. Detection antibody: 82943-4-PBS. Standard: Eg0147. Range: 0.313-10 ng/mL



Cytometric bead array standard curve of MP00006-2, PDGFRB Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82943-5-PBS. Detection antibody: 82943-4-PBS. Standard: Eg0147. Range: 1.25-20 ng/mL