

For Research Use Only

DGAT1 Recombinant antibody, PBS Only



Catalog Number: 82945-1-PBS

Basic Information

Catalog Number: 82945-1-PBS	GenBank Accession Number: BC015762	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 8694	CloneNo.: 230091E1
Source: Rabbit	UNIPROT ID: O75907	
Isotype: IgG	Full Name: diacylglycerol O-acyltransferase homolog 1 (mouse)	
Immunogen Catalog Number: AG32717	Calculated MW: 488 aa, 55 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 82945-1-RR (DGAT1 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 82945-1-PBS in a different storage buffer formulation.



WB result of DGAT1 antibody (82945-1-RR; 1:3000; incubated at room temperature for 1.5 hours) with sh-Control and sh-DGAT1 transfected HeLa cells. This data was developed using the same antibody clone with 82945-1-PBS in a different storage buffer formulation.