

For Research Use Only

IFT140 Recombinant antibody, PBS Only (Detector)



Catalog Number: 82947-5-PBS

Basic Information

Catalog Number: 82947-5-PBS	GenBank Accession Number: BC035577	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 9742	CloneNo.: 230231A12
Source: Rabbit	UNIPROT ID: Q96RY7	
Isotype: IgG	Full Name: intraflagellar transport 140 homolog (Chlamydomonas)	
Immunogen Catalog Number: AG11119	Calculated MW: 1462 aa, 165 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00029-1, IFT140 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82947-2-PBS. Detection antibody: 82947-5-PBS. Standard: Ag11119. Range: 2.5-20 ng/mL.