

For Research Use Only

HEBP1 Recombinant antibody, PBS Only (Capture)



Catalog Number: 82950-5-PBS

Basic Information

Catalog Number: 82950-5-PBS	GenBank Accession Number: BC016277	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 50865	CloneNo.: 230156F10
Source: Rabbit	UNIPROT ID: Q9NRV9	Affinity: $K_D = 4.79 \times 10^{-10}$
Isotype: IgG	Full Name: heme binding protein 1	$K_{off} = 3.22 \times 10^{-4}$
Immunogen Catalog Number: AG10538	Calculated MW: 189 aa, 21 kDa	$K_{on} = 6.72 \times 10^5$

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00032-2, HEBP1 Recombinant Matched Antibody Pair - PBS only. Capture antibody: 82950-5-PBS. Detection antibody: 82950-4-PBS. Standard: Ag10538. Range: 1.25-20 ng/mL