

For Research Use Only

# EFNA2 Recombinant antibody, PBS Only (Detector)



Catalog Number:82951-4-PBS

## Basic Information

<b>Catalog Number:</b> 82951-4-PBS	<b>GenBank Accession Number:</b> NM_001405	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 1mg/ml	<b>GeneID (NCBI):</b> 1943	<b>CloneNo.:</b> 230251D6
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> O43921	
<b>Isotype:</b> IgG	<b>Full Name:</b> ephrin-A2	
	<b>Calculated MW:</b> 213 aa, 24 kDa	

## Applications

**Tested Applications:**  
Indirect ELISA, Cytometric bead array

**Species Specificity:**  
Human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00081-1, EFNA2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82951-5-PBS. Detection antibody: 82951-4-PBS. Standard: SY00429. Range: 0.313-10 ng/mL