

For Research Use Only

uPAR, PLAUR Recombinant antibody, PBS Only (Detector)



Catalog Number: 82953-3-PBS

Basic Information

Catalog Number:

82953-3-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG28056

GenBank Accession Number:

BC002788

GeneID (NCBI):

5329

UNIPROT ID:

Q03405

Full Name:

PLAUR

Calculated MW:

37 kDa

Purification Method:

Protein A purification

CloneNo.:

230148E10

Applications

Tested Applications:

Cytometric bead array, ELISA

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Standard curve of MP00129-1 measured by sandwich ELISA. Capture antibody: 82953-2-PBS. Detection antibody: 82953-3-PBS. Range: 125-8000 pg/mL