

For Research Use Only

# uPAR, PLAUR Recombinant antibody, PBS Only (Capture)



Catalog Number: 82953-4-PBS

## Basic Information

Catalog Number:

82953-4-PBS

Size:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG28056

GenBank Accession Number:

BC002788

GeneID (NCBI):

5329

UNIPROT ID:

Q03405

Full Name:

PLAUR

Calculated MW:

37 kDa

Purification Method:

Protein A purification

CloneNo.:

230148F8

## Applications

Tested Applications:

Cytometric bead array, ELISA

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Standard curve of MP00129-2 measured by sandwich ELISA. Capture antibody: 82953-4-PBS. Detection antibody: 82953-2-PBS. Range: 1.25-8000 pg/mL



Standard curve of MP00129-2 measured by cytometric bead array. Capture antibody: 82953-4-PBS. Detection antibody: 82953-2-PBS. Range: 1.25-20 ng/mL