

For Research Use Only

MAX Recombinant antibody, PBS Only (Capture)

Catalog Number: 82958-4-PBS



Basic Information

| | | |
|--|--|---|
| Catalog Number: 82958-4-PBS | GenBank Accession Number: BC003525 | Purification Method: Protein A purification |
| Size: 1mg/ml | GeneID (NCBI): 4149 | CloneNo.: 230215C10 |
| Source: Rabbit | UNIPROT ID: P61244 | Affinity: $K_D=4.96 \times 10^{-10}$ |
| Isotype: IgG | Full Name: MYC associated factor X | $K_{Off}=1.93 \times 10^{-4}$ |
| Immunogen Catalog Number: AG0680 | Calculated MW: 18 kDa | $K_{On}=3.90 \times 10^5$ |

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00126-2, MAX Recombinant Matched Antibody Pair - PBS only. Capture antibody: 82958-4-PBS. Detection antibody: 82958-3-PBS. Standard: Ag0680. Range: 0.625-20 ng/mL