

For Research Use Only

SCN3B Recombinant antibody, PBS Only (Detector)

Catalog Number: 82959-5-PBS



Basic Information

Catalog Number:

82959-5-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG23201

GenBank Accession Number:

BC126265

GeneID (NCBI):

55800

UNIPROT ID:

Q9NY72

Full Name:

sodium channel, voltage-gated, type

III, beta

Purification Method:

Protein A purification

CloneNo.:

230316H5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00132-2, SCN3B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82959-2-PBS. Detection antibody: 82959-5-PBS. Standard: Ag23201. Range: 0.625-20 ng/mL