

For Research Use Only

SCN3B Recombinant antibody, PBS Only (Detector)

Catalog Number: 82959-6-PBS



Basic Information

Catalog Number: 82959-6-PBS	GenBank Accession Number: BC126265	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 55800	CloneNo.: 230316H6
Source: Rabbit	UNIPROT ID: Q9NY72	
Isotype: IgG	Full Name: sodium channel, voltage-gated, type III, beta	
Immunogen Catalog Number: AG23201		

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00132-1, SCN3B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82959-2-PBS. Detection antibody: 82959-6-PBS. Standard: Ag23201. Range: 0.625-20 ng/mL