

For Research Use Only

SCN3B Recombinant antibody, PBS Only (Detector)

Catalog Number: 82959-8-PBS



Basic Information

Catalog Number: 82959-8-PBS	GenBank Accession Number: BC126265	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 55800	CloneNo.: 230316B7
Source: Rabbit	UNIPROT ID: Q9NY72	
Isotype: IgG	Full Name: sodium channel, voltage-gated, type III, beta	
Immunogen Catalog Number: AG23201		

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP00132-4, Human SCN3B Recombinant Matched Antibody Pair - PBS only. 82959-2-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag23201. 82959-8-PBS was HRP conjugated as the detection antibody. Range: 0.156-10 ng/mL



The mean SCN3B concentration was determined to be 32.36 ng/mL in HEK-293 cell extract based on a 3.40 mg/mL extract load.