

For Research Use Only

RAB32 Recombinant antibody, PBS Only (Detector)

Catalog Number: 82960-2-PBS



Basic Information

Catalog Number: 82960-2-PBS	GenBank Accession Number: BC015061	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 10981	CloneNo.: 230315B11
Source: Rabbit	UNIPROT ID: Q13637	
Isotype: IgG	Full Name: RAB32, member RAS oncogene family	
Immunogen Catalog Number: AG1450	Calculated MW: 25 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00043-1, RAB32 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82960-4-PBS. Detection antibody: 82960-2-PBS. Standard: Ag1450. Range: 2.5-20 ng/mL.