

For Research Use Only

ARHGEF19 Recombinant antibody, PBS Only (Capture)



Catalog Number: 82966-3-PBS

Basic Information

Catalog Number: 82966-3-PBS	GenBank Accession Number: BC040640	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 128272	CloneNo.: 230290B5
Source: Rabbit	UNIPROT ID: Q8IW93	
Isotype: IgG	Full Name: Rho guanine nucleotide exchange factor (GEF) 19	
Immunogen Catalog Number: AG11645	Calculated MW: 802 aa, 89 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00072-1, ARHGEF19 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82966-3-PBS. Detection antibody: 82966-4-PBS. Standard: Ag11645. Range: 0.625-20 ng/mL.



Cytometric bead array standard curve of MP00072-2, ARHGEF19 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82966-3-PBS. Detection antibody: 82966-5-PBS. Standard: Ag11645. Range: 0.625-20 ng/mL.