

For Research Use Only

ACAP1 Recombinant antibody, PBS Only (Detector)

Catalog Number: 82967-3-PBS



Basic Information

Catalog Number: 82967-3-PBS	GenBank Accession Number: BC018543	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 9744	CloneNo.: 230370G9
Source: Rabbit	UNIPROT ID: Q15027	
Isotype: IgG	Full Name: ArfGAP with coiled-coil, ankyrin repeat and PH domains 1	
Immunogen Catalog Number: AG34309	Calculated MW: 740 aa, 82 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP01212-1, ACAP1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82967-2-PBS. Detection antibody: 82967-3-PBS. Standard: Ag34309. Range: 1.25-20 ng/mL