

For Research Use Only

DDX60 Recombinant antibody, PBS Only (Capture)



Catalog Number: 82969-2-PBS

Basic Information

Catalog Number: 82969-2-PBS	GenBank Accession Number: BC038115	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 55601	CloneNo.: 230293B9
Source: Rabbit	UNIPROT ID: Q8IY21	
Isotype: IgG	Full Name: DEAD (Asp-Glu-Ala-Asp) box polypeptide 60	
Immunogen Catalog Number: AG34406	Calculated MW: 1712 aa, 198 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00038-1, DDX60 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82969-2-PBS. Detection antibody: 82969-3-PBS. Standard: Ag34406. Range: 1.25-20 ng/mL.