

For Research Use Only

# CTGF Recombinant antibody, PBS Only (Detector)



Catalog Number: 82976-2-PBS

## Basic Information

Catalog Number:

82976-2-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM\_001901

GeneID (NCBI):

1490

UNIPROT ID:

P29279

Full Name:

connective tissue growth factor

Calculated MW:

38 kDa

Purification Method:

Protein A purification

CloneNo.:

230324C3

## Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

Human

## Background Information

### Storage

Storage:

Store at -80°C.

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Sandwich ELISA standard curve of MP00091-1, CTGF  
Recombinant Matched Antibody Pair, PBS Only.  
Capture antibody: 82976-4-PBS. Detection  
antibody: 82976-2-PBS. Standard: SY01606. Range:  
0.313-10 ng/mL