

For Research Use Only

LIMK2 Recombinant antibody

Catalog Number: 82977-1-RR



Basic Information

Catalog Number: 82977-1-RR	GenBank Accession Number: BC013051	Purification Method: Protein A purification
Size: 1000 µg/ml	GeneID (NCBI): 3985	CloneNo.: 230283D10
Source: Rabbit	UNIPROT ID: P53671	Recommended Dilutions: IF 1:125-1:500
Isotype: IgG	Full Name: LIM domain kinase 2	
Immunogen Catalog Number: AG3011	Calculated MW: 686 aa, 78 kDa	

Applications

Tested Applications: IF/ICC, ELISA	Positive Controls: IF : HepG2 cells,
Species Specificity: Human	

Background Information

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed HepG2 cells using LIMK2 antibody (82977-1-RR, Clone: 230283D10) at dilution of 1:250 and CoraLite®488-Conjugated AffiniPure Goat Anti-Mouse IgG(H+L) (SA00013-1).