

For Research Use Only

# POU4F3 Recombinant antibody, PBS Only (Detector)



Catalog Number: 82980-2-PBS

## Basic Information

<b>Catalog Number:</b> 82980-2-PBS	<b>GenBank Accession Number:</b> BC104923	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 1mg/ml	<b>GeneID (NCBI):</b> 5459	<b>CloneNo.:</b> 230345B9
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> Q15319	
<b>Isotype:</b> IgG	<b>Full Name:</b> POU class 4 homeobox 3	
<b>Immunogen Catalog Number:</b> AG16004	<b>Calculated MW:</b> 338 aa, 37 kDa	

## Applications

**Tested Applications:**  
Indirect ELISA, Cytometric bead array

**Species Specificity:**  
Human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Cytometric bead array standard curve of MP00087-1, POU4F3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82980-4-PBS. Detection antibody: 82980-2-PBS. Standard: Ag16004. Range: 0.625-20 ng/mL.



Cytometric bead array standard curve of MP00087-2, POU4F3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82980-5-PBS. Detection antibody: 82980-2-PBS. Standard: Ag16004. Range: 0.625-20 ng/mL.