

For Research Use Only

POU4F3 Recombinant antibody, PBS Only (Capture)



Catalog Number: 82980-4-PBS

Basic Information

Catalog Number: 82980-4-PBS	GenBank Accession Number: BC104923	Purification Method: Protein A purification
Size: 1mg/ml	GeneID (NCBI): 5459	CloneNo.: 230345D12
Source: Rabbit	UNIPROT ID: Q15319	
Isotype: IgG	Full Name: POU class 4 homeobox 3	
Immunogen Catalog Number: AG16004	Calculated MW: 338 aa, 37 kDa	

Applications

Tested Applications:
Indirect ELISA, Cytometric bead array

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00087-1, POU4F3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82980-4-PBS. Detection antibody: 82980-2-PBS. Standard: Ag16004. Range: 0.625-20 ng/mL.



Cytometric bead array standard curve of MP00087-3, POU4F3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82980-4-PBS. Detection antibody: 82980-3-PBS. Standard: Ag16004. Range: 0.625-20 ng/mL.