

For Research Use Only

OGDHL Recombinant antibody, PBS Only (Capture)



Catalog Number: 82981-5-PBS

Basic Information

Catalog Number: 82981-5-PBS	GenBank Accession Number: BC026320	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 55753	CloneNo.: 230228H6
Source: Rabbit	UNIPROT ID: Q9ULDO	
Isotype: IgG	Full Name: oxoglutarate dehydrogenase-like	
Immunogen Catalog Number: AG10794	Calculated MW: 1010 aa, 115 kDa	

Applications

Tested Applications:
Cytometric bead array, ELISA

Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.

Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Standard curve of MP00023-2 measured by cytometric bead array. Capture antibody: 82981-5-PBS. Detection antibody: 82981-4-PBS. Range: 2.5-20 ng/mL