

For Research Use Only

MMP3 Recombinant antibody, PBS Only (Capture)



Catalog Number: 82984-1-PBS

Basic Information

Catalog Number:

82984-1-PBS

Size:

1mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC074815

GeneID (NCBI):

4314

UNIPROT ID:

P08254

Full Name:

matrix metalloproteinase 3
(stromelysin 1, progelatinase)

Purification Method:

Protein A purification

CloneNo.:

230535B3

Applications

Tested Applications:

Indirect ELISA, Cytometric bead array

Species Specificity:

Human

Background Information

Storage

Storage:

Store at -80°C.

Storage Buffer:

100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve ofMP00083-1, MMP3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82984-1-PBS. Detection antibody: 82984-2-PBS. Standard: Eg0382. Range: 0.625-20 ng/mL



Cytometric bead array standard curve ofMP00083-3, MMP3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82984-1-PBS. Detection antibody: 82984-4-PBS. Standard: Eg0382. Range: 0.625-20 ng/mL