

For Research Use Only

ADAMDEC1 Recombinant antibody, PBS Only (Capture)

Catalog Number: 82985-2-PBS



Basic Information

| | | |
|---|--|---|
| Catalog Number: 82985-2-PBS | GenBank Accession Number: BC074877 | Purification Method: Protein A purification |
| Size: 1mg/ml | GeneID (NCBI): 27299 | CloneNo.: 230367B10 |
| Source: Rabbit | UNIPROT ID: O15204 | |
| Isotype: IgG | Full Name: ADAM-like, decysin 1 | |
| Immunogen Catalog Number: AG12265 | Calculated MW: 470 aa, 53 kDa | |

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP00007-3, ADAMDEC1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82985-2-PBS. Detection antibody: 82985-5-PBS. Standard: Ag12265. Range: 0.625-20 ng/mL