

For Research Use Only

ADAMDEC1 Recombinant antibody, PBS Only (Detector)



Catalog Number: 82985-3-PBS

Basic Information

Catalog Number: 82985-3-PBS	GenBank Accession Number: BC074877	Purification Method: Protein A purification
Size: 1 mg/ml	GeneID (NCBI): 27299	CloneNo.: 230367B11
Source: Rabbit	UNIPROT ID: O15204	
Isotype: IgG	Full Name: ADAM-like, decysin 1	
Immunogen Catalog Number: AG12265	Calculated MW: 470 aa, 53 kDa	

Applications

Tested Applications:
Indirect ELISA
Species Specificity:
Human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage Buffer:
100% PBS pH 7.3

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data