

For Research Use Only

# Mouse Hif1a Recombinant antibody, PBS Only (Capture)



Catalog Number: 82989-1-PBS

## Basic Information

<b>Catalog Number:</b> 82989-1-PBS	<b>GenBank Accession Number:</b> NM_001313919	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 1mg/ml	<b>GeneID (NCBI):</b> 15251	<b>CloneNo.:</b> 230101A3
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> Q61221	
<b>Isotype:</b> IgG	<b>Full Name:</b> hypoxia inducible factor 1, alpha subunit	
<b>Immunogen Catalog Number:</b> AG33620	<b>Calculated MW:</b> 93 kDa	

## Applications

**Tested Applications:**  
Sandwich ELISA, Indirect ELISA

**Species Specificity:**  
Mouse

## Background Information

### Storage

**Storage:**  
Store at -80°C.

**Storage Buffer:**  
100% PBS pH 7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Sandwich ELISA standard curve of MP00056-2, Mouse H1fa Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 82989-1-PBS. Detection antibody: 82989-3-PBS. Standard: Ag33620. Range: 0.156-10 ng/mL.